



MentorMed
IMPLANT



www.mentormedimplant.com

**MADE IN
GERMANY**

Modern Medical Equipment LLC

Founded in 2006, Headquartered in Dubai, U.A.E, **MENTORMED** is showing the dental implant industry how a strong business model can service customers effectively and economically. **MENTORMED** is a full line top quality, compatible dental implants and accessories at a reasonable cost. we supply Dental Clinics, Doctors, Hospitals and laboratories with every type of specialized dental equipment for endodontics, implantology, Bone and membrane, and dental kits. we built a well-known, impeccable reputation in the community. Our management has many years of experience in the field of dental implants.


At **MENTORMED** we are not only the most competitively priced, service-oriented supply company in the nation; we are now the best full-line compatible dental implants and accessories supply company. We are located in Dubai but serve customers across the country of United Arab Emirates and beyond.

Two key differences between **MENTORMED** and most other distributors are our broad product offering and our supply base. We are very proud to represent the leaders in dental implant manufacturing, and we work closely with our suppliers to make sure that we continue to meet the needs of our customers.

True to our entrepreneurial spirit, we continue to seek new solutions and products that meet the new demands of the sector. Therefore, we entered into an agreement with ALPHA DENT IMPLANTS GMBH a German manufacturer for the exclusive production and manufacture of **MENTORMED** made of premium high-strength materials of high quality and bio compatible. as our motto goes "where quality meets service", our promise to our customers is to continue improvement and sourcing out the best quality products available, at reasonable cost combined with the best service that goes beyond the sale.



TESTED BENEFITS



1 Platform Switching

- Switching Platform
- Narrow connection of the implant prevents compression in cortical layer

2 Micro thread

- Micro threads for better contact with a cortical plate

3 Intermediate thread

- Bone thread on the corticospongoid level

4 Active thread

- Has multiple levels of height, which allows enlarging the surface of the implant
- Allows achieving reliable initial stability in soft bone

5 Antirotational sulcus

6 3D Active surface

- Unique active and hydrophilic surface of a new generation

7 Apical part

- Curving in apical part prevents the damage of Sinus

MENTORMED IMPLANT




CONVENTIONAL SIGNS

Reference		Implant diameter	
Implant Name	30 MM	Implant length	
Expiry date	medical device REF XXXXXX	Lot number	
Barcode	Implant XXXX Dm=XXmm; Ln=XXmm;		
Expiry date	XXXX-XX LOT XXXXXXXXXXXX		
Lot number	(01) XXXXXXXXXXXXXXXX		
Serial number	(17) XXXX-XX-XX (11) XXXX-XX-XX		
Manufacturer address	(10) XXXXXXXXXXXX		
	(21) XXXXXXXXXXXXXXXXXXXXXXXX		
	Alpha Dent Implants GmbH Hanauer Str.1-5, 75181 Pforzheim, Germany		Certification in compliance with Appendix II, clause 3 of Guideline 93/42 EWG 93/42 EWG Medical Products Council
	Made in Germany		Website
	www.alpha-dent.net		
			

STERILITY STERILE R


All production of MENTORMED Implants is thoroughly sterilized. Double sterile packing ensures sterility maintenance. Due to a high technical level of production and the strictest control of sterility, MENTORMED Implants meets the requirements of the highest quality standards, confirmed by appropriate certificates.

Positioning of the implant in bone




Bone depth on the side of cortical plate must be not less than 1 mm.

Distance between the implant and tooth root



The distance between an implant and a tooth root must be minimum 1.5 mm, i.e. the distance from the implant site center to the next tooth must be about the $\frac{1}{2}$ diameter of the implant's body + 1.5 mm.



Distance depends on diameter of implants installed



The distance between two implants must be a sum of $\frac{1}{2}$ diameters of the body of the implants installed + 3 mm (distance between implants).

SLA SURFACE

SLA surface represents the latest generation of dental implant surface treatment technology. This surface is active and hydrophilic; it has distinct multi-layered microporous structure. The Thickness of oxide layer of the implant surface reaches 10–12 microns. The entire thickness of the oxide layer is imbued with micro pores that combine with each other and create a multi-layered surface. This way, the bone tissue grows deeper both in and between micro pores. Thereby, implant surface area is significantly increased and, thus the contact between bone tissue and implant surface is strengthened. Titanium oxide film is enriched with calcium hydroxyapatite, and the implant surface is contaminant-free. It has prominent hydrophilic characteristics (absorption).




Researches show that SLA surface stimulates active bone growth over the implant surface that results in mechanical stability. This surface property allows performing orthopedic rehabilitation at an earlier stage and significantly reduces the risk when the implant is installed under the one-stage procedure.

AVAILABLE DIAMETERS AND LENGTHS

PLATFORM 3.3

Ø Size (mm)	Length (mm)						
	4.0	6.75	8	10	11.5	13	16
3.5							
4.0							
4.5							
5.5							



Implant has a platform with an inner cone and a hexagon for positioning. This type of platform allows creating a conical connection with suprastructures. This connection is more stable and airtight. 3.25 mm diameter implants are installed in front and lateral parts (canine teeth, incisors, premolar teeth).

1- Length of the subcortical part of the implant is 4.0 mm.
Full implant length is 6.2 mm.

Ø Size (mm)	Length (mm)	Reference	
3.5	10	IMM 3.5x10	
	11.5	IMM 3.5x11.5	
	13	IMM 3.5x13	
	16	IMM 3.5x16	
4.0	4.0	IMM 4.0x4.0	
	6.75	IMM 4.0x6.75	
	8.0	IMM 4.0x8.0	
	10	IMM 4.0x10	
	11.5	IMM 4.0x11.5	
	13	IMM 4.0x13	
	16	IMM 4.0x16	
4.5	4.0	IMM 4.5x4.0	
	6.75	IMM 4.5x6.75	
	8.0	IMM 4.5x8.0	
	10	IMM 4.5x10	
	11.5	IMM 4.5x11.5	
	13	IMM 4.5x13	
5.5	4.0	IMM 5.5x4.0	
	6.75	IMM 5.5x6.75	
	8.0	IMM 5.5x8.0	
	10	IMM 5.5x10	
	11.5	IMM 5.5x11.5	
	13	IMM 5.5x13	

1– Length of the subcortical part of the implant is 4.0 mm.
Full implant length is 6.2 mm.



Cover screw*



Platf.	Length (L).	Reference
	4.70	GLYRT


Anodised








*Screw included with each implant.








Metric



Table of sequential use of conical diamond like coated (DLC) drills for implant Mentormed





Drill diameter, mm	1.5	2.0	2.8	3.0	3.45	3.9	4.7	6.0
rate, RPM	D1 800-1200 D2 800-1200 D3 800-1000 D4 800-1000	D1 800-1000 D2 800-1000 D3 700-900 D4 600-800	D1 600-800 D2 600-800 D3 500-600 D4 400-500	D1 500-700 D2 600-800 D3 500-600 D4 400-500	D1 500-700 D2 600-800 D3 500-600 D4 400-500	D1 400-600 D2 500-700 D3 500-600 D4 400-500	D1 400-600 D2 400-700 D3 500-600 D4 400-500	D1 400-600 D2 400-700 D3 500-600 D4 400-500

	Implants 3.5mm - bone type I, II
3.5	
	Implants 3.5mm - bone type III, IV
3.5	
	Implants 4.0mm - bone type I, II
4.0	
	Implants 4.0mm - bone type III
4.0	
	Implants 4.0mm - bone type IV
4.0	

Diamond like coated drills in Implants system are developed with the aim to form bone cavity for implants Mentormed because the drill shape recreates perfectly the shape of the implant without thread.

**Table of sequential use of conical diamond
like coated (DLC) drills for implant Mentormed**

Drill diameter, mm	1.5	2.0	2.8	3.0	3.45	3.9	4.7	6.0
rate, RPM	D1 800-1200 D2 800-1200 D3 800-1000 D4 800-1000	D1 800-1000 D2 800-1000 D3 700-900 D4 600-800	D1 600-800 D2 600-800 D3 500-600 D4 400-500	D1 500-700 D2 600-800 D3 500-600 D4 400-500	D1 500-700 D2 600-800 D3 500-600 D4 400-500	D1 400-600 D2 500-700 D3 500-600 D4 400-500	D1 400-600 D2 400-700 D3 500-600 D4 400-500	D1 400-600 D2 400-700 D3 500-600 D4 400-500

	Implants 4.5mm - bone type I, II
4.5	
	Implants 4.5mm - bone type III, IV
4.5	
	Implants 5.5mm - bone type I, II, III
5.5	
	Implants 5.5mm - bone type IV
5.5	

Diamond like coated drills in Implants system are developed with the aim to form bone cavity for implants Mentormed because the drill shape recreates perfectly the shape of the implant without thread.


ABUTMENTS

DIRECT-TO-IMPLANT RECONSTRUCTIONS



SYMBOLOLOGY

 Rotatory element	 Tx30 connection	 Made from steel
 Non-rotatory element	 Size in millimetres	 Made from cobalt-chromium + castable plastic
 Use with manual torque	 45° screw support	 Made from cobalt-chromium
 Maximum operating torque	 90° screw support	 Made from PEEK
 Ratchet torque range	 Use in rotation with a CA	 Made from plastic
 Galaxy connection	 Maximum rotation speed	 Recommended sterilisation temperature
 Screw connection	 Maximum number of uses	 Unsterilised product
 Kirator connection	 Single-use product	 Use with abundant irrigation
 Nature connection	 Made from grade 5 ELI (extra-low interstitial) titanium	 Maximum angle
 Basic connection	 Made from grade 2 titanium	
 XDrive connection	 Made from stainless steel	



Ø 10




Code **PMAKM3.0**




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Pre - Milled Abutment Konus 3



Ø 6.0




Code **FFAM 3x2** **FFAM 3x5**




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Free Form Abutment



Ø 3.35



Code **TAKM 3 7** **TAKM 3 9** **TAKM 3 12**




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Abutment Konus 3-7mm, 9mm, 12mm

Straight titanium antirotational conical abutment. Abutment includes a screw.
Available in three sizes: 7 mm, 9 mm, and 12 mm.



Ø 5.0



Code **TAWKM 3 9**




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Abutment Wide 3 9mm

Wide antirotational titanium conical abutment with hexagon. Available sizes 9 mm. Mainly used in the area of posterior teeth.





NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Abutment Konus 3-1mm, 2mm, 3mm

Anatomical titanium shoulder abutment straight conical. Available in three shoulder sizes: 1.5, 2.5, 3.5 mm.

Indicator includes 0.5 mm of subcortical immersion.




NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Abutment Wide Konus 3-1mm, 2mm, 3mm

Anatomical titanium shoulder wide abutment straight conical. Available in three shoulder sizes: 1.5, 2.5, 3.5 mm.

Indicator includes 0.5 mm of subcortical immersion.




NARROW CONICAL PLATFORM 2.1 mm (3.0)

Angulated Titanium Abutment 15deg 3 Konus-1mm, 2mm, 3mm

Angular anatomical titanium conical abutment 15° with shoulder 1.5–3.5 mm.

It is used for aesthetic tasks solving in vestibular area, and replicating the shape of gingival line.


Indicator includes 0.5 mm of subcortical immersion.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Angulated Titanium Abutment Konus 3 25 deg






NARROW CONICAL PLATFORM 2.1 mm (3.0)

Angulated Titanium Abutment 25 deg 3 Konus-1mm, 2mm, 3mm


Indicator includes 0.5 mm of subcortical immersion.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Angulated Titanium Abutment Konus 3 15 deg

Angular antirotational titanium conical abutment 15° without a shoulder. It is used in cases, when the implant is installed at an angle. This type of the abutment allows creating the aesthetically correct construction with less effort for technician and orthopedist.




NARROW CONICAL PLATFORM 2.1 mm (3.0)

Ball Attachment Konus 3

D Ø 4.2

d Ø 2.5




Code BAKM 3/2 BAKM 3/3 BAKM 3/4 BAKM 3/5



Metal Cap

Silicone Cap




NARROW CONICAL PLATFORM 2.1 mm (3.0)

Locator 3

D Ø 1.25

d Ø 3.9



Code LCKM3/2 LCKM3/3 LCKM3/4 LCKM3/5 LCKM3/6



Metal Cap




A structure with conical ball attachment is used for fixation of removable denture. A cap made of stainless steel can be included, as well as an intermediate silicone cap of three different degrees of hardness: soft, standard, and strong. Ball attachment set does not include caps. It can be purchased separately.




silicone caps

Plastic disc



Ø 4.5



Code

PTAKM3




Screw for the abutment


NARROW CONICAL PLATFORM 2.1 mm (3.0)

Composed Hex Abutment Konus 3

Titanium-based plastic conical abutment is used with single implants.



Ø 4.2



Code


PAKM3 NH




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Plastic Cylinder Abutment Konus 3 Non Hex



Ø 3,4



Code


TBKM3




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Base Konus 3



Ø 1.5



Code

TBKM3-SC




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Base Konus 3






NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Platform for CAD/CAM Konus Non Hex 3

D Ø 4.6


d Ø 3.3



Code **TPCKM3-NH** **TPCKM3-NH-1** **TPCKM3/NH-2** **TPCKM3/NH-3**



Retaining screw



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Locator 3

D Ø 4.6

NON HEX


d Ø 3.3



Code **TPCKM3** **TPCKM3/1** **TPCKM3/2** **TPCKM3/3**



Retaining screw



Code




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)



NARROW CONICAL PLATFORM 2.1 mm (3.0)



Code



Screw for the abutment




HEALING CAPS FOR IMPLANTS



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Base Konus 3




Code	TBKM3
------	-------

Narrow MentorMed Healing Cap provides perfect gingival margin. It fits for any following options of dental prosthesis.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Healing Cap Konus 3



Code	HCKM 3/2	HCKM 3/3	HCKM 3/4	HCKM 3/5	HCKM 3/7
------	----------	----------	----------	----------	----------


Titanium conical healing cap is installed on the implant and is intended to form gingival line and gingival papilla. It is available in following sizes: 2, 3, 4, 5, 7 mm. The size is defined according to the height of gingiva.

Indicator includes 0.5 mm of subcortical immersion.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Wide Healing CAP




Code	HCWKM 3/2	HCWKM 3/3	HCWKM 3/4	HCWKM 3/5	HCWKM 3/7
------	-----------	-----------	-----------	-----------	-----------

Indicator includes 0.5 mm of subcortical immersion.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Wide Healing CAP




Code	HCWKM 6.5/3	HCWKM 6.5/4	HCWKM 6.5/5	HCWKM 6.5/7
------	-------------	-------------	-------------	-------------

Indicator includes 0.5 mm of subcortical immersion.



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Wide Healing CAP



Code	HCWKM 8.0/3	HCWKM 8.0/4	HCWKM 8.0/5	HCWKM 8.0/7
------	-------------	-------------	-------------	-------------

Indicator includes 0.5 mm of subcortical immersion.



TRANSFERS AND ANALOGS FOR **IMPLANTS**



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Implant Analog Konus 3




ANKM3

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Digital Implant Analog Konus 3



D-ANKM3



Ø 3,5




Code




Screw for the abutment

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Scan Body Implant CAD/CAM PEEK



Ø 3,5



Code


ATNM9




Retaining screw

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Scan Abutment for CAD/CAM Conical



Ø 3,5



Code

ATNM9-Ti



Retaining screw


NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium Abutment Training Narow 9mm



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer Konus Open Tray 3 13mm




Code	TCKM3/13
------	----------



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer Konus Open Tray 3 15mm



Code	TCKM3/15
------	----------





Retaining screw



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer Konus Close Tray 3 9mm



Code	TCKM3/9
------	---------




Retaining screw



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer Konus Open Tray 3 13mm



Code	TCKM3/13
------	----------





Retaining screw



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer Konus Open Tray 3 15mm



Code	TCKM3/15
------	----------



Retaining screw





The length of conical closed tray transfer is 13.75 mm. It is used in closed tray impression coping.

MULTI UNITS FOR IMPLANTS





NARROW CONICAL PLATFORM 2.1 mm (3.0)
ONE -PIECE ANGULAR MULTI-UNIT 17°



Code	MUBKM3/1701	MUBKM3/1702	MUBKM3/1703
------	-------------	-------------	-------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)
ONE -PIECE ANGULAR MULTI-UNIT 30°
Indicator includes 0.5 mm of subcortical immersion.




Code	MUBKM3/3001	MUBKM3/3002	MUBKM3/3003
------	-------------	-------------	-------------



Multi-unit system is used for a prosthesis with screw retention. One-piece angular multi-unit allows rehabilitation at non-parallel implants by correcting the angle. Angular multi-unit is available in options 17° and 30° options

One-piece angular multi-unit 45° conical. Multi-unit system is used for a prosthesis with screw retention. One-piece angular multi-unit allows rehabilitation at non-parallel implants by correcting the angle. The new angular multi-unit is available in 45° option.



Ø 4.3





Code	
------	--




NARROW CONICAL PLATFORM 2.1 mm (3.0)
ONE -PIECE ANGULAR MULTI-UNIT 45°

Ø 4.3



Code	
------	--






Code **AEAKM3/001+** **AEAKM3/002+** **AEAKM3/003+**

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Antirotational Esthetic Abutment Konus 3

Aesthetic abutment (straight multi-unit) with a plastic sleeve is used in order to restore a single structure or a bridgewith a screw retention. Screw in with a maximum force of 35N



D Ø 3.3



d Ø 4.0

Code **ANAEMM**

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Analog for angular multi-unit 17°/30°


It serves as an analog of an implant with an **aesthetic antirotation abutment (straight multi-unit)** in the dental technical lab. Allows preparing laboratory models.



Code **D-ANAEMM**

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Digital Analog for Straight Multi-Unit



Ø 4.5



Code **TCAEMM**



Retaining


NARROW CONICAL PLATFORM 2.1 mm (3.0)

Transfer for angular multi-unit 17°/30°t



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Base for Out-Burnt Abutment




Code	CCMM Base
------	-----------



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Titanium platform for CAD_CAM MU





Code	TPCM-MU
------	---------



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Healing Cap PEEK Konus 3



Code	HCPKM3/2.25
------	-------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Healing Cap for Angular Multi-Unit 17°/30°



Code	HCAEMM
------	--------

Titanium healing cap for multi-unit is intended to form gingival line and gingival papilla



NARROW CONICAL PLATFORM 2.1 mm (3.0)

Healing Cap wide for Angular Multi-Unit 17°/30°



Code	HCAEMM- W
------	-----------

Titanium healing cap for multi-unit is intended to form gingival line and gingival papilla



NARROW CONICAL PLATFORM 2.1 mm (3.0)
Plastic cuff multi-unit 17°/30°



Code	PMMM
------	-------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)
Temporary PEEK cuff multi-unit 17°/30°



Code	PEMMM
------	--------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)
Titanium cuff multi-unit 17°/30°



Code	TSAAMMT
------	----------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)
Scan Abutment for MU



Code	ATM-MU
------	---------------



NARROW CONICAL PLATFORM 2.1 mm (3.0)
Temporary abutment PEEK
straight multi unit



PEMKM

NARROW CONICAL PLATFORM 2.1 mm (3.0)
Plastic cuff straight multi-unit 17°/30°



PMKM

NARROW CONICAL PLATFORM 2.1 mm (3.0)
Titanium cuff multi-unit 17°/30°



TSAAKM


NARROW CONICAL PLATFORM 2.1 mm (3.0)
Healing Cap for Aesthetic Abutment Konus



HCAEKM

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Antirotational Esthetic Abutment Konus 3





4.7
D Ø
2.1



Code	AEAKM3/001	AEAKM3/002	AEAKM3/003
------	------------	------------	------------



Retaining screw




Code	TCAEKM
------	--------

NARROW CONICAL PLATFORM 2.1 mm (3.0)


Open Tray Transfer for Aesthetic Antirotational Abutment

Open tray transfer for aesthetic antirotational abutment (straight multi-unit) with a size of 10 mm. It is used in open tray impression coping.



D Ø 3.3

d Ø 4.0



Code	ANAEKM
------	--------

NARROW CONICAL PLATFORM 2.1 mm (3.0)

Analog of Aesthetic Antirotation Abutment (Straight Multi-Unit)

It is used as an implant analog with multi-unit in the dental laboratory. It allows preparing laboratory models.



SURGICAL KIT



DRILL

Diamond like coated drills

Cylindrical



TCAM 1.5



TCAM 2.0



TCAM 2.8

Tapered finishing



TCAM 3.0



TCAM 3.45



TCAM 3.9



TCAM 4.7



TCAM 5.7

COUNTERSINKS



TDCSM3.2



TDCSM3.75



TDCSM4.2



TDCSM5



TDCSM6

Drill extension



TDM-E

RR

Torque ratchet



RRM



STOPPERS

for 2.8 drills and 2.0

6 mm | 8 mm | 10 mm | 11.5 mm | 13 mm



for TCA 3.0, 3.45, 3.9, 4.7, 5.7 drills

6 mm | 8 mm | 10 mm | 11.5 mm | 13 mm



IMPLANT DRIVERS



HPIDKM 3/9



HPIDKM 3/18



IDKM3/9



IDKM3/18

Universal wrench for abutment screw



HDSM9



HDSM18

Direction indicator depth gauge



PPM



PPM 2.8

SURGICAL GUIDE KIT



DRILL

Diamond like coated drills



GDM-1.9/6.75 GDM-1.9/8 GDM-1.9/10 GDM-1.9/11.5 GDM-1.9/13 GDM-1.9/16



GDM-2.8/6.75 GDM-2.8/8 GDM-2.8/10 GDM-2.8/11.5 GDM-2.8/13 GDM-2.8/16



GDM-3.2/6.75 GDM-3.2/8 GDM-3.2/10 GDM-3.2/11.5 GDM-3.2/13 GDM-3.2/16



GDM-3.6/6.75 GDM-3.6/8 GDM-3.6/10 GDM-3.6/11.5 GDM-3.6/13 GDM-3.6/16



GDM-4.0/6.75 GDM-4.0/8 GDM-4.0/10 GDM-4.0/11.5 GDM-4.0/13 GDM-4.0/16



GDM-4.8/6.75 GDM-4.8/8 GDM-4.8/10 GDM-4.8/11.5 GDM-4.8/13 GDM-4.8/16

MOUNTER PLATFORM

RATCHET
ADAPTER

HANDPIECE
ADAPTER

DRILL FOR
FIXATION PIN

DRILL
STOPPERS



MPM Konus 3+ MPM Konus 3+2 MPM Konus 3+4 MPM Konus 3+6



ADLM Hex 4 ADSM Hex 4



MDM Hex.3



PDM-1.5



DSM - H1 DSM - H1.5

GUIDE
SLEEVES

EXTRACTORS

FIXATION PIN

TISSUE
PUNCH

PILOT DRILL
BONE PROFILE

DRILL
EXTENSION

UNIVERSAL
SCREWDRIVER

IMPLANT
DRIVER



STM-1.5

GSM-D5-h4.0

SEM 1.72

PTM-1.5

G5/TP5.0

GBMM 4.8

TDM-E

HDSM9

HDSM18

HPIDKM 3/18

IMPLANT GUIDES FOR SYSTEMS

Straumann | Neodent | Nobel | Megagen | NeoBiotech | Osstem | DIO | Bauer's

Suitable for most of the most popular implant systems

Different types of offset: +, +2, +4, +6



N/cm ²	Product name	Possible length				Compatibility with systems
50 N	Mounter Platform Type I	+	+2	+4	+6	Megagen AnyOne, NeoBiotech, Osstem (4.0+), DIO (4.0+), Bauer's
50 N	Mounter Platform Type II	+	+2	+4	+6	Megagen AnyRidge
40 N	Mounter Platform Type III	+	+2	+4	+6	Nobel CC NP
50 N	Mounter Platform Type IV	+	+2	+4	+6	Nobel CC RP
40 N	Mounter Platform Type V	+	+2	+4	+6	Straumann BLT NC
50 N	Mounter Platform Type VI	+	+2	+4	+6	Straumann BLT RC
40 N	Mounter Platform Type VII	+	+2	+4	+6	Megagen Mini
50 N	Mounter Platform Type VIII	+	+2	+4	+6	Neodent
40 N	Mounter Platform Type IX	+	+2	+4	+6	Nobel Replace NP
50 N	Mounter Platform Type X	+	+2	+4	+6	Straumann BLX



World wide Distributor

**MODERN
MEDICAL
EQUIPMENT**

Modern Medical Equipment L.L.C

Tel: +971 4 236 6604

AL Khaleej Road, Dnaat Khaleej Bldg.

Shop # 15, P.O. Box: 251395 Dubai, UAE

dubai@mmeco.ae

www.mmeco.ae



Manufacturer By
Alpha Dent Implants GmbH
Pforzheim
Hanauer Str. 8, 75181
Germany



Visit our website



www.mentormedimplant.com



info@mentormedimplant.com

**CE
1023**